


DETERGENT - DISINFECTANT  
MEDICAL DEVICE  
CLASS IIa



# UMONIUM<sup>38</sup>

## MEDICAL TISSUES



✓ FORMULA FREE FROM  
→ CMR COMPONENT\*\*  
→ ENDOCRINE DISRUPTOR

“ Instant disinfection to hand ”

SPORICIDE



**UMONIUM<sup>38</sup> MEDICAL TISSUES** is the first cleaning and sporicide disinfecting wipe. It kills germs within 5 seconds of contact time. It can be used to clean and disinfect fragile materials such as fibre optic, polycarbonate, latex, vaginal probes, endoscopes, unit external operating unit instruments, etc. Ultra-resistant and moist, each wipe disinfects an area of 3 m<sup>2</sup>.

\*\* Carcinogenic, Mutagenic and Reprotoxic





# UMONIUM<sup>38</sup> MEDICAL TISSUES

1 wipe (20 x 20 cm) = 3 m<sup>2</sup>



✓ **FORMULA FREE FROM**  
→ CMR COMPONENT\*\*  
→ ENDOCRINE DISRUPTOR

## INDICATIONS

Cleaning and disinfection for non-invasive medical devices.

## COMPOSITION

Woven (box of 100 and individual pouches) and non-woven (box of 200) wipes soaked in a solution of **UMONIUM<sup>38</sup> EQUIPMENTS**: N-benzyl-N-dodecyl-N, N-dimethyl-ammonium chloride/N-benzyl-N, N-diméthyl-N-tetradecyl-ammonium chloride. 11,9g/L.

## REGULATORY STANDARDS

- Bactericidal, fungicidal and sporicidal activity of the soaked wipe is validated according to standards AFNOR NF T72-190, ASTM 2967-15, EN 16615
- The bactericidal, fungicidal, levuricidal, mycobactericidal and virucidal activity of the impregnation solution is validated according to EN 13727, EN 14561, EN 13624, EN 14562, EN 14348, EN 14563 and EN 14476 standards
- Complies with ED 93/42/EEC relating to medical devices (annexe II, excl. §4)
- Quality and traceability guaranteed by ISO 9001:2015/EN ISO 13485:2016

## INSTRUCTIONS FOR USE

**Cleaning/Pre-disinfection:** Take out a wet wipe. Wet wipe to eliminate all dirt. Rinse thoroughly with water\* or using a cloth dampened with water\*. **Disinfection:** Take out a new wet wipe. Wipe the clean devices for at least 5 seconds. Leave to dry or rinse with water\* (or using a cloth dampened with water\*) if necessary. Make sure to fully close the cover after use.

## PROPERTIES

- Medical device Class IIa
- Detergent - disinfectant wipe with a broad spectrum of activity
- Non-toxic, no dangerous vapours
- LD<sub>50</sub> (rats): > 5840 mg/kg
- Non-corrosive: pH 7
- Ready-to-use product not irritating to the skin (no transdermal toxicity)
- Aldehyde, peroxide and biguanide free, no other carcinogenic, mutagenic, reprotoxic components
- Non-coagulant
- Non-flammable
- Bactericide, fungicide and sporicidal activities
- 100% compatible on fibre optic, rubber, polycarbonate, acrylic compound, glass, Pyrex, PU, PVC, HDPE, PET, neoprene, latex, silicone, paints and stainless steel 410
- Effective, fast and efficient (dissolves blood)
- Effective in the presence of organic materials
- Works within 5 seconds
- Remains on unrinse surfaces (bacteriostatic film)
- No residue after rinsing
- Stability: 36 months date of manufacture
- Biodegradable solution and recyclable wipe



ITEM No.	NAME	BOX
PF 10760	Box of 100 wipes (20 x 20 cm)	6 x 100
PF 10780	Refill of 100 wipes for dispensing box	12 x 100
PF 10762	Box of 200 wipes (13 x 14.5 cm)	6 x 200
EMB 5206	Stainless steel bracket	1 item
PF 10748	Box of 10 individual pouches	36 x 10
PF 10761	500 individual pouches	500 x 1

\* Clean drinking water at room temperature  
\*\* Carcinogenic, Mutagenic and Reprotoxic

[www.huckerts.net](http://www.huckerts.net)

FP00052 EN - Ed. 3.3 - 02/2019



**BELGIUM**  
20 Avenue Lavoisier  
1300 Wavre  
Tel.: +32 (0)67 89 41 00  
Fax: +32 (0)67 84 37 67

**FRANCE**  
38 Avenue du Centre  
78230 Le Pecq-sur-Seine  
Tel.: +33 (0)1 3976 1505  
Fax: +33 (0)1 3976 1115

**G.D. LUXEMBOURG**  
19 Rue de l'Industrie  
8069 Bertrange  
Tel.: +352 26 39 42 60  
Fax: +352 26 39 42 70


**ITALY**  
20 Via Giulio Natta  
05100 Terni  
Tel.: +39 0744 80 03 87  
Fax: +39 0744 81 73 39

**SINGAPORE**  
51 Changi Business Park Central 2  
The Signature, #04-05  
Singapore 486066  
Tel.: +65 983 45 333



UMONIUM<sup>38</sup> is a registered trademark of Laboratoire Huckert's International

## PROCEDURE FOR USE



### STEP 1 : CLEANING/PRE-DISINFECTION



Take out a wet wipe

Wet wipe to eliminate all dirt


Rinse thoroughly with water\* or using a cloth dampened with water\*

### STEP 2 : DISINFECTION



Take out a new wet wipe

Wipe the clean devices for at least 5 seconds



Leave to dry **OR** rinse with water\* (or using a cloth dampened with water\*) if necessary